Create it



### SURG

BONE SURGERY P05-P12
SCRAPER P13-P14
EXTRACTION P15-P16
SOCKET LIFT (CRESTAL APPROACH) P17-P20
IMPLANT PREPARATION P21-P24
SOCKET LIFT (CRESTAL APPROACH) P25-P26
SINUS LIFT P27
SINUS MEMBRANE DETACHMENT P28

#### PERIO

SCALING	P29-P31
PERIO	P32
PERIO-CONTROL	P33-P34

#### ENDO

ENDODONTICS	P35-P36
TIP KIT	P37-P40

\*Pictures provided by Yasuhiro Kizu, D.D.S., Ph.D. OMFC (ORAL & MAXILLOFACIAL CARE CLINIC YOKOHAMA)

## **VARIOSURG TIP LINEUP**

SURG mode for bone surgery

PERIO mode for maintenance

ENDO mode for root canal treatment



The tips are designed to preserve tissue, maintain cutting speed and easily cut into bone. Saw like teeth on an uneven blade edge enhances bone cutting. This also reduces clogging to maintain cutting speed.

**BONE SURGERY** 



IMPLANT PREPARATION



PERIO



SCRAPER



SINUS LIFT



PERIO-CONTROL



**EXTRACTION** 



SINUS MEMBRANE DETACHMENT



**ENDODONTICS** 



SOCKET LIFT (CRESTAL APPROACH)



SCALING



TIP KIT







## VarioSurg 3

## **EVOLVED POTENTIAL**

#### Embodying the Pursuit of Perfection

Ultrasonic bone surgery systems are essential in modern surgical and implant treatment.

The powerful VarioSurg3 handpiece features a slim and lightweight body
that ensures minimal hand fatigue and outstanding visibility. A variety of tips makes procedures easier
and stress-free to help provide more precise and efficient treatment.





## **BONE SURGERY**

30.4

MODEL H-SG1

POWER LEVEL SURG 150%

SG1A

POWER LEVEL SURG 80%

For use in deeper regions requiring longer tips.

31.0



ORDER CODE

Z305101

- TiN coating Dots mark 3, 6 and 9 mm from tip end
- Five teeth: 0.6 mm thick

HIGH POWER

31.0



ORDER CODE

Z305151

- TiN coating Dots mark 3, 6 and 9 mm from tip end
- Five teeth: 0.6 mm thick

33.6



ORDER CODE Z305138



• TiN coating • Dots mark 3, 6, 9, 12 and 15 mm from tip end



MODEL SG8

30.4

POWER LEVEL SURG 80%

ORDER CODE

Z305109

MODEL H-SG8

HIGH POWER


POWER LEVEL SURG 150%

ORDER CODE

Z305155

MODEL SG8A

For use in deeper regions requiring longer tips.



• TiN coating • Dots mark 3, 6, 9, 12 and 15 mm from tip end

• TiN coating • Dots mark 3, 6 and 9 mm from tip end

• Three teeth: 0.6 mm thick

• Three teeth: 0.6 mm thick

• TiN coating • Dots mark 3, 6 and 9 mm from tip end

• Three teeth: 0.6 mm thick



ORDER CODE Z305122



• TiN coating • 90 degree right curved tip\*

Five teeth; 0.6 mm thick

• Five teeth: 0.6 mm thick



Anterior mandib ular harvesting SG14R, SG14L





ORDER CODE

Z305123

- TiN coating 90 degree left curved tip\*



HIGH POWER





ORDER CODE

Z305156



• Three teeth; 0.6 mm thick

MODEL H-SG8R POWERLEVEL SURG 150% MODEL H-SG8L

POWER LEVEL SURG 150%







ORDER CODE

Z305157

• TiN coating • 90 degree left curved tip\*

• Three teeth; 0.6 mm thick

\* The direction of the tip's curve is defined by the tip's anterior view.

• TiN coating • Right curved tip\*

• Three teeth: 0.6 mm thick

30.0

ORDERCODE


Z305102




ORDER CODE

Z305103

- TiN coating Left curved tip\*
- Three teeth: 0.6 mm thick







ORDER CODE Z305135

31.0

• TiN coating • Dot mark 3 mm from tip end

• Five teeth: 0.8 mm thick



**SG30** 

corticotomies.

For narrow cutting to reduce bone loss. Also suitable for

27.8

POWER LEVEL SURG 50%

MODEL SG58

POWER LEVEL SURG 80%

For removing tori on the posterior mandible or mandibular lingual surface.





ORDER CODE

Z305137

• TiN coating • Sharp edge • 0.5 mm thick



ORDER CODE

Z305141

• TiN coating • Dots mark 3, 6 and 9 mm from tip end

Three teeth: 0.6 mm thick



The cutting part is bone torus

\* The direction of the tip's curve is defined by the tip's anterior view.

12

SCRAPER

#### For flattening the alveolar crest or harvesting bone fragments.

ORDER CODE

Z305105

**SCRAPER** 

MODEL

SG3

POWER LEVEL SURG 80%

MODEL SG4

POWER LEVEL SURG 80%


MODEL SG5

POWER LEVEL SURG 80%

ORDER CODE

Z305106

For removal of diseased granulation tissues from extraction sites.













ORDER CODE

Z305104

• TiN coating • Triple sided edge spatula tip Single-edged scraper



Double-edged scraper

• TiN coating • Edge spatula tip

• TiN coating • Rounded edge spatula tip

Single-edged scraper







Ruptures the periodontal ligament to streamline suture removal. Also effective in extracting ankylosed teeth.

**EXTRACTION** 

**SG17** 

• TiN coating • 0.7 mm thick

POWER LEVEL SURG 80%

MODEL SG18R

POWER LEVEL SURG 50%

MODEL SG18L

POWER LEVEL SURG 50%

ORDER CODE







Z305132

ORDER CODE





ORDER CODE

Z305133

Z305134


27.7

• TiN coating • Left angled tip\* • 0.7 mm thick



• TiN coating • Right angled tip\* • 0.7 mm thick






# SOCKET LIFT (CRESTAL APPROACH)

## SOCKET LIFT (CRESTAL APPROACH) Socket lift procedures.

The socket lift tip enhances the safety of procedures, being less invasive and enabling sinus lifts with short chair times.

Osteotomy with ultrasonic vibration reduces the risk of damage to the Schneiderian membrane.


Eliminating hammering taps through the osteotome technique reduces patient discomfort and potential dizziness.



A starting point is added to trigger bone formation.

Bone formation to prepare implant site.

Implant site expansion and formation.



Maxillary sinus perforation and bone formation on alveolar ridge.

Reliably elevate the maxillary sinus membrane with internal irrigation.

Forms implant preparation site once maxillary sinus membrane is elevated.

19 20

IMPLANT PREPARATION For preparation of implant sites, adding starting points, and other steps.

For preparation of implant sites, flattening bone surfaces, and other steps. IMPLANT PREPARATION

MODEL SG15A

POWER LEVEL SURG 50%

SG15B

POWER LEVEL SURG 50%





28.5



ORDER CODE

Z305124

Z305125

ORDER CODE

- Diamond coating Tip-end diameter; 0.7 mm
- Markings every 5 mm and 10 mm from tip end

- Diamond coating Tip-end diameter; 1.3 mm
- Markings every 5 mm and 10 mm from tip end



Implant preparation SG15A, SG15B

SG15C

POWER LEVEL SURG 50%

MODEL SG15D

POWER LEVEL SURG 50%







ORDER CODE

Z305126

- Diamond coating Tip-end diameter; 0.9 mm
- Markings every 5 mm, 10 mm and 12 mm from tip end





ORDER CODE

Z305127

- Diamond coating Tip-end diameter; 1.3 mm
- Marking every 5 mm, 10 mm and 12 mm from tip end



Implant preparation SG15C, SG15D



22

IMPLANT PREPARATION For expanding and forming implant preparation site when approaching from the alveolar crest.

SG16A

POWER LEVEL SURG 50%

SG16B

POWER LEVEL SURG 50%





ORDER CODE

Z305128

- · Diamond coating · Diameter; 2.0 mm
- Markings every 5 mm and 10 mm from tip end



29.9



ORDER CODE

Z305129

- Diamond coating Diameter; 2.6 mm
- Markings every 5 mm and 10 mm from tip end



Implant preparation

SCL2D

POWER LEVEL SURG 50%

MODEL SCL4D

POWER LEVEL SURG 80%

INTERNAL IRRIGATION





ORDER CODE Z305182

INTERNAL IRRIGATION





ORDER CODE Z305184

Diamond coating



· Diamond coating







Formation of implant preparation site SCL2D, SCL4D

ORDER CODE

Z305170

Before the mucous membrane of maxillary sinus is elevated SCL1, SCL3, SCL5

ORDER CODE

Z305172

INTERNAL IRRIGATION

29.8

POWER LEVEL SURG 50%

MODEL SCL3 MODEL SCL5 POWER LEVEL SURG 50% POWER LEVEL SURG 50%

INTERNAL IRRIGATION

SCL3D

SOCKET LIFT (CRESTAL APPROACH)

For procedures such as osteoloniy of maximary so crest and implant preparation site adjustments.

Diamond coating



INTERNAL IRRIGATION 28.1 ORDER CODE

POWER LEVEL SURG 80%

Z305175

SCL5D

· Diamond coating

For procedures such as osteotomy of maxillary sinus approaching from the alveolar

POWER LEVEL SURG 80%

· Diamond coating

using markings on the tip.

MODEL SCL1D

INTERNAL IRRIGATION

28.7

depth can be measured

ORDER CODE

Z305171

POWER LEVEL SURG 80%

Cutting the cavity floor of implant preparation site SCL1D, SCL3D, SCL5D

Implant preparation site depth can be measured using markings on the tip.



MODEL SCL1

INTERNAL IRRIGATION

ø2.4 <

28.7





MODEL SG6D

MODEL SG9

POWER LEVEL SURG 50% MODEL SG10

POWER LEVEL SURG 50% MODEL SG11

POWER LEVEL SURG 50%



29.2



ORDER CODE

POWER LEVEL SURG 80%

Z305107

ORDER CODE

POWER LEVEL SURG 50%

Z305108

Diamond coating

MODEL SG7D



• Flat circular convex elevator • Angled at 90°

24.0

Z305110

• Flat circular convex elevator • Angled at 135°





ORDER CODE

Z305111

ORDER CODE

Z305112

Cone compressor

31.4



· Diamond coating

Avulsion of wisdom teeth and impacted teeth




Lateral maxillary posterior bone augmentation

ORDER CODE

Sinus membrane detachment SG9, SG10



Sinus membrane elevation



## SCALING

MODEL G1-S

POWER LEVEL SURG 50%

- 0.8 G1-S



ORDERCODE

Z305113

Single rooted / Multi rooted

MODEL P20-S

POWER LEVEL PERIO 100%

MODEL P25R-S

POWER LEVEL PERIO 100%

MODEL P25L-S

POWER LEVEL PERIO 100%

Removal of Biofilm furcation area.

0.8 P20-S



36.4

ORDER CODE

Z305114





ORDER CODE

Z305115

Removal of Biofilm furcation area.






ORDER CODE

Z305116

Left curved type\*



Right curved type\*



Single rooted / Multi rooted P25R-S, P25L-S

\*The direction of the tip's curve is defined by the tip's anterior view.

POWER LEVEL PERIO 100%

MODEL V-P10



For the maintenance of prostheses and implants. Produced from special material to allow safe debridement without causing damage to metals and implants.









ORDER CODE

Z305117

ORDER CODE

Y900184

- Includes E-Tip replacement wrench
- Plastic Tip is not included

- Pack of 3
- V10-S holder is not included



**V-P11R** 

or contours.

MODEL V-P11L

Double bent tips. This tip is extra long to easily access difficult to reach areas and to avoid interference with neighbouring teeth

A longer version of the popular V-P10 tip.

MODEL V-P12















ORDER CODE

Y1002165

ORDER CODE Y1002166 ORDER CODE Y1002167

- Pack of 3 Right curved type\*
- V10-S holder is not included

- Pack of 3 Left curved type\* V10-S holder is not included

 Pack of 3 V10-S holder is not included





33

<sup>■</sup> V-P11R/L, V-P12 can be used only for VarioSurg 3 \*The direction of the tip's curve is defined by the tip's anterior view.

• Diamond coating


• For posterior teeth (right angled\*)

**ENDODONTICS** 

MODEL E30RD-S POWERLEVEL ENDO 50% MODEL E30LD-S POWERLEVEL ENDO 50% MODEL E31D-S POWERLEVEL ENDO 50%













ORDER CODE

Z305119



ORDER CODE

Z305120



ORDER CODE

Z305118

Diamond coating

For posterior teeth (left angled\*)

· Diamond coating

• For anterior and posterior teeth (70°)


Diamond coating

• For anterior teeth (90°)



## ORDERCODE Z305121

<sup>\*</sup> The direction of the tip's angle is defined by the tip's anterior view.



TIP KIT

### MODEL Basic H-S Kit

#### Contents

- H-SG1, SG3, SG5, SG6D, SG7D, SG11
- Tip holder



SINUS MEMBRANE DETACHMENT

ORDER CODE

Y1002775

### MODEL Bone Cut Kit

#### Contents

- SG1, SG2R, SG4, SG2L, SG6D
- Tip holder



ORDER CODE

Y900688

#### BONE SURGERY

- SG15A, SG15B, SG16A, SG16B
  - Tip holder

Contents



MODEL Implant Preparation Kit

ORDER CODE

IMPLANT PREPARATION

Y900774

### MODEL Socket Lift (CRESTAL APPROACH) Kit for Regular Size Implant

#### Contents

- SCL1D, SCL2D, SCL1
- VS Tip Wrench Tip holder Manual



ORDER CODE

SOCKET LIFT (CRESTAL APPROACH

Y1002841

### MODEL Sinus Lift Kit

#### Contents

- SG1, SG3, SG6D, SG9, SG10, SG11
- Tip holder



ORDERCODE Y900689

### MODEL Endo-S Kit

#### Contents

- G1-S, E30RD-S, E30LD-S, E31D-S, E32D-S
- Tip holder



ORDER CODE

Y900691

• VS Tip Wrench • Tip holder • Manual

• SCL2D, SCL3D, SCL4D, SCL3

Contents



MODEL Socket Lift (CRESTAL APPROACH) Kit for Wide Size Implant

SOCKET LIFT (CRESTAL APPROACH)

Y1002842

ORDER CODE

## VarioSurg 3

### Advantages

- Super slim handpiece with LED
- Extensive tip range for a wide range of applications
- Link 2 system with Surgic Pro2 (implant motor)





#### Complete Set

•	VarioSurg3 (230V)	Y1002726
Optic	MODEL	ORDER CODE

- Control Unit LED Handpiece with 2m cord FC-78 Foot Control

- Irrigation tube (5 pcs.) and other accessories
   Basic H-S Kit (H-SG1, SG3, SG5, SG6D, SG7D, SG11 and Tip holder)

#### Complete Set without Foot Control

Optic	VarioSurg3 (230V) Non FT	Y1002248

- Control Unit LED Handpiece with 2m cord Sterilization Cassette
- Handpiece Stand Irrigation tube (5 pcs.) and other accessories
   Basic H-S Kit (H-SG1, SG3, SG5, SG6D, SG7D, SG11 and Tip holder)

#### Specifications (Control Unit)

- Frequency: 28-32 kHz
- Power Supply : 230V 50/60 Hz
- Irrigation Flow Rate : 10-75 mL / min
- Programs : SURG x 5, ENDO x 2, PERIO x 2
- Dimensions : W 265 x D 220 x H 103 mm

NAKANISHI INC. www.nsk-dental.com 700 Shimohinata, Kanuma, Tochigi 322-8666, Japan NSK Oceania Pty. Ltd. www.australia.nsk-dental.com Unit 12, 809-821 Botany Road, Rosebery NSW 2018, Australia

NSK Oceania Ltd. www.new-zealand.nsk-dental.com Ground Floor, Suite D, 272 Parnell Road, Parnell, Auckland 1052, New Zealand

Specifications are subject to change without notice.