

# THE NEW STYLE ELECTRIC THE COMPACT MAKES A BIG CHANGE

## NLZ E ELECTRIC MICROMOTOR SYSTEM



# NEW PRODUCT

A totally new installation style only from NSK



## Specification comparison

		Our existing model	NLZ • NLZ E
Motor specs	Number of rotations	2,000-40,000 min <sup>-1</sup>	<b>100-40,000 min<sup>-1</sup></b>
	Max torque	3.4 Ncm	<b>4.2 Ncm</b> <span style="border: 1px solid orange; border-radius: 50%; padding: 2px;">120% High Torque</span>
Corresponding model to endodontic mode		-	<b>○ (NLZ E)</b>

## NLZ-E

### Electric Micromotor

PORTABLE SYSTEM WITH ENDO FUNCTION

### + FREE NLZ Endo Handpiece and ISO-E Lux Probe

- Easily convert air driven handpieces to electric with the all new NLZ-E.
- Correctly transmit reciprocating actions for a smooth endodontic treatment.
- New NLZ ENDO 6:1 reciprocating handpiece.
- Endo torque calibration function – This function measures the 'no-load' current values and corrects the torque accordingly, resulting in greatly enhance safety in endodontic treatment.
- Contra-Check Function helps detect abnormalities caused by deterioration of the contra-angle handpiece and insufficient maintenance before use.

MODEL	ORDER CODE	RRP PRICE
NLZ-E Portable	Y1500361	<b>\$5,501</b>
NLZ Endo Handpiece	C1130	<b>FREE</b> RRP-\$827
ISO E-Lux Probe	Z1265	<b>FREE</b> RRP-\$688



## NLZ E





## NLX nano

Electric Micromotors

PORTABLE SYSTEM WITH NANO CONTROL

This system can help any user of air driven handpieces to convert to an electric system.



MODEL	ORDER CODE	PROMO PRICE
NLX Portable	Y1001362	<b>\$4,112</b> RRP \$5,140



## INTEGRATED SYSTEM WITH MULTIPAD CONTROL

This system can be easily integrated into your existing chair to control up to 2 electric micromotors and 1 Varios scaling system.

MODEL	ORDER CODE	PROMO PRICE
NLX nano standard with MultiPad	NLXNANOIntegratedpad	<b>\$3,507</b> RRP \$4,383
NLX plus with MultiPad	NLXPLUSIntegratedpad	<b>\$3,691</b> RRP \$4,614



## Varios 170 iPiezo engine



## Varios 170

LED Illumination

INTEGRATED MULTI FUNCTION ULTRASONIC SCALER

The smallest module for easy integration into new and existing dental equipment. Contents: Built in module, Handpiece, Handpiece cord, 3 scaling tips (G4, G6, G8), 3 x Tip replacement wrench with torque limiter, Sterilisation case. Installation not included.

MODEL	ORDER CODE	PROMO PRICE
Varios 170 Optic Complete Set	Y1001351	<b>\$1,672</b> RRP \$2,090
Varios 170 Non-Optic Complete Set	Y1001352	<b>\$1,187</b> RRP \$1,484



20% OFF

# S-Max M Handpieces

Advanced Stainless Steel



MODEL **S-Max M900L**  
ORDER CODE P1254  
• 26W Power  
• Standard Head  
**PROMO PRICE**  
**\$1,023**  
RRP \$1,279



MODEL **S-Max M800L**  
ORDER CODE P1253  
• 23W Power  
• Miniature Head  
**PROMO PRICE**  
**\$1,023**  
RRP \$1,279



MODEL **S-Max M65**  
ORDER CODE H1008  
• Non-Optic  
• 1:1 Direct Drive  
• Single spray  
• HP/CA burs  
**PROMO PRICE**  
**\$605**  
RRP \$757

## Mini Air Micro Motors



MODEL **S-Max M205**  
ORDER CODE M1007001  
• Non-Optic  
• 1:1 Direct Drive  
• Internal Spray  
• Non retraction valve  
**PROMO PRICE**  
**\$899**  
RRP \$1,124

Non-optic options available. Other connections offered for KaVo, Sirona, W&H, Bien Air.



MODEL **S-Max M95L**  
ORDER CODE C1023  
• Optic  
• Increasing 1:5  
• Quattro spray  
• FG Burs  
**PROMO PRICE**  
**\$1,306**  
RRP \$1,632



MODEL **S-Max M25L**  
ORDER CODE C1024  
• Optic  
• 1:1 Direct Drive  
• Single spray  
• CA burs  
**PROMO PRICE**  
**\$1,002**  
RRP \$1,252



MODEL **S-Max M15L**  
ORDER CODE C1025  
• Optic  
• 4:1 Reduction  
• Single Spray  
• CA Burs  
**PROMO PRICE**  
**\$1,220**  
RRP \$1,525

"The Osseo 100 unit has been extremely easy and affordable to use. Being able to compare RFA readings gives me confidence that we are ready to progress to the restorative phase."



Dr David Attia BOH (DentSci) (Griffith), GradDipDent (Griffith), MSc (Oral Implantology) (Goethe), PGDipClin (Orth)(CoL), FICOI

## Osseo 100 + FREE Sterile covers

Osseo 100 measures implant stability and osseointegration using RFA, to enhance decisions about when to load the implant.

MODEL **Osseo 100**  
ORDER CODE Y1004175  
**RRP PRICE**  
**\$3,025**

MODEL **Osseo 100 Sterile Covers**  
ORDER CODE 55105  
**PROMO PRICE**  
**FREE**  
RRP \$99

FREE Sterile Covers THREE PACKS OF 20



## Simple and precise root canal treatment

### Endodontic Motors

#### ENDO-MATE TC2

MODEL **ENDO-MATE TC2 (230V)**  
**MP-F16R head**  
ORDER CODE Y1001027

- Large LCD display with 5 programs for different file systems
- Supports most major brands of Ni-Ti files
- Torque control facility ranging from 0.1 to 4.0 NCM depending on head:gear ratio

**PROMO PRICE**  
**\$1,676**  
RRP \$2,095



#### ENDO-MATE AT

MODEL **Endo-Mate AT (230V)**  
ORDER CODE Y1002190

- Fits all major brands of Files
- Torque control, with integrated auto-reverse function
- Compact stylish control/monitor pad for simple and user friendly operation

**PROMO PRICE**  
**\$1,340**  
RRP \$1,675



### Apex Locator

#### IPEX II

MODEL **iPEX II Complete set**  
ORDER CODE Y1002208

- Accurate measurement of any root canal
- SmartLogic technology
- Operates in all conditions, including dry, wet and bloody canals

**PROMO PRICE**  
**\$1,256**  
RRP \$1,570



20% OFF

"The NSK Endomate TC2 has been my motor of choice for the past 10 years. It is a well balanced, cordless motor that delivers consistent torque and speed for the endodontic treatment."



Dr Morgan Wu  
Specialist Endodontist  
Bds (Otago), Dcd (Endo)(Melb)

# Varios Combi Pro Basic Kit

## + FREE Perio Kit

### 2 Maintenance Systems in 1

Varios Combi Pro combines ultrasonic scaling and powder therapy for a more effective and comfortable oral hygiene and biofilm control

Varios Combi Pro

#### Basic Kit

MODEL **Varios Combi Pro Basic Kit**  
ORDER CODE Y1002843

- Advanced optic iPiezo scaling for General, Perio & Endo treatment.
- Powder Therapy to clean Supragingival areas efficiently and quickly.
- User friendly handpiece and nozzle for easy access.

#### Perio Kit

MODEL **Varios Combi Pro Perio Kit**  
ORDER CODE Y1003042

- Remove biofilm from Subgingival areas with the gentle glycine based Perio-Mate Powder.
- Sterile single-use plastic nozzle tips.
- Safely remove biofilm inside periodontal pockets and around dental implants.



**FREE PERIO KIT WORTH \$2,628**

MODEL	ORDER CODE	RRP PRICE
Varios Combi Pro Basic Kit	Y1004194	<b>\$7,274</b>
Varios Combi Pro Perio Kit	Y1003042	<b>FREE</b> RRP-\$2,628



Scan the QR code for more information

## Varios2 Scaler Handpiece

MODEL	ORDER CODE	PROMO PRICE
Varios2 Lux HP Optic Scaler Handpiece	E351050	<b>\$1,370</b> RRP-\$1,644
Varios2 HP Non-Optic Scaler Handpiece	E350050	<b>\$834</b> RRP-\$982



**15% OFF VARIOS HANDPIECES**

## Wide choice of Ultrasonic Scaler Tips

### Buy 4 tips + FREE torque wrench

MODEL	RRP PRICE
Ultrasonic Scaler Tips	From <b>\$128</b>

MODEL	ORDER CODE	PROMO PRICE
Varios Torque Wrench	Z221076	<b>FREE</b> RRP-\$122



Scan the QR code for full tip guide



**BUY 4 TIPS GET FREE TORQUE WRENCH**

## VIVA ace + FREE ISO E-Luxe Probe

### MOBILE DENTISTRY SYSTEM

For more comfortable and more advanced home visits for dental treatment. VIVA ace combines the functionality of a dental treatment centre in a compact body with multipurpose micromotor, ultrasonic scaler and more. Its lightweight build means it is easy to carry and straight forward to set up in different locations. The first of its kind with high-functionality, supports more comfortable dental treatment wherever the patient.

Price excludes handpieces.

MODEL	ORDER CODE	PROMO PRICE
VIVA ace complete set	Y1004228	<b>\$17,550</b> RRP \$19,500



## ISO E-Lux Probe

### INTENSE DIRECT ILLUMINATION

Intense direct illumination for diagnosis of cracks in teeth plus occlusal and interproximal caries not visible to the naked eye.

MODEL	ORDER CODE	PROMO PRICE
ISO E-Lux Probe	Z1265	<b>\$550</b> RRP \$688

**FREE**  
with  
VIVA ace  
purchase

VIVA ace  
MOBILE DENTISTRY SYSTEM

Ask about  
our bundle  
options



# AUTOCLAVE

## iClave plus

### B-TYPE VACUUM AUTOCLAVE

- 18 litre chamber capacity.
- 20% more sterilisation space than a conventional stainless steel chamber of the same size.
- Copper chamber minimises internal temperature fluctuations.
- Super fast cycle as quick as 18 minutes, including drying.
- One year warranty.

Unit size: H428 x W445 x D532mm, chamber 18L.

MODEL	ORDER CODE	PROMO PRICE
iClave plus	Y1003077	<b>\$8,910</b> RRP \$9,900

**10%  
OFF**

iClave plus



NSK Oceania Pty Ltd [www.nskoceania.com.au](http://www.nskoceania.com.au)

Product Information | Technical Support

Freephone: 1300 44 33 21 (Australia) Email: [info@nskoceania.com.au](mailto:info@nskoceania.com.au)