

# Special Offers

July - September 2023



## Powerful perfection

Maximum power 150%\*

Stabilising the power balance for more efficient procedures



## VarioSurg 3

Ultrasonic Bone Surgery System

\*Compared to our existing model

# HANDPIECES

## High Speed Handpieces

Ti-Max Series and S-Max

BUY 3 GET 1 FREE\*

### Ti-Max Z

Titanium Contra Angle

- Optic
- 1:5 Speed Increasing
- Clean Head System (Anti-Retraction)
- Titanium body with DURAGRIP
- Cellular Glass Optics
- Push Button Chuck
- Microfilter
- 2 Years warranty



### Switch Spray Function

- Switch between Mist Spray and Jet Spray
- Mist Spray (Water and air)
- Jet Spray (Water only)
- No equipment modification required

MODEL	Ti-Max Z95L
ORDER CODE	C1034SW
RRP	\$2,121
SAVE	\$2,121



### Ti-Max Z

Titanium Turbine

Power 26W

- Optic
- Clean Head System (Anti-Retraction)
- Torque: 26W (model Z900L)
- Smallest head : ø12.5 x H 13.1 mm
- Speed: 320,000~400,000 min<sup>-1</sup>
- Titanium body with DURAGRIP
- Quattro Water Spray
- Push Button chuck
- Ceramic bearings
- Microfilter
- 2 Years warranty



MODEL	Ti-Max Z900L
ORDER CODE	P1111
RRP	\$1,981
SAVE	\$1,981

MODEL	Ti-Max Z900KL
ORDER CODE	P1113
RRP	\$2,074
SAVE	\$2,074



MODEL	Ti-Max Z900SL
ORDER CODE	P1115
RRP	\$2,074
SAVE	\$2,074

MODEL	Ti-Max Z900WL
ORDER CODE	P1121
RRP	\$2,074
SAVE	\$2,074

### S-Max M

Stainless Steel Turbine

Power 26W

- Optic
- Clean Head System (Anti-Retraction)
- Torque: 26W (model M900L)
- Speed: 320,000~400,000 min<sup>-1</sup>
- Stainless Steel Body
- Quattro Water Spray
- Push Button chuck
- Ceramic bearings
- Microfilter
- Simple Cartridge Replacement
- 1 Year warranty

MODEL	S-Max M900L
ORDER CODE	P1254
RRP	\$1,381
SAVE	\$1,381

MODEL	S-Max M900KL
ORDER CODE	P1258
RRP	\$1,591
SAVE	\$1,591



MODEL	S-Max M900SL
ORDER CODE	P1262
RRP	\$1,591
SAVE	\$1,591

MODEL	S-Max M900WL
ORDER CODE	P1266
RRP	\$1,591
SAVE	\$1,591

\*Buy 3 Get 1 free offer available on listed handpieces only. Free handpiece given on the equal or lesser value handpiece. Mix and match accepted on the listed handpiece.

## S-MAX 20% OFF

MODEL	S-Max M25L
ORDER CODE	C1024

- Optic
- 1:1 Direct Drive
- Single Spray
- For CA burs

RRP	\$1,353
SALE PRICE	\$1,082
SAVE	\$271



MODEL	S-Max M25
ORDER CODE	C1027

- Non-Optic
- 1:1 Direct Drive
- Single Spray
- For CA burs

RRP	\$949
SALE PRICE	\$759
SAVE	\$190



MODEL	Ti-Max Z25L
ORDER CODE	C1038

- Optic
- 1:1 Direct Drive
- Single Spray
- For CA burs (ø2.35)

RRP	\$1,491
SALE PRICE	\$1,193
SAVE	\$298



## NSK LED COUPLING 50% OFF

## Enhance your visability with NSK LED Couplings



MODEL	PTL-CL-LED
ORDER CODE	P1001600

- Optic
- Stainless Steel Body
- Titanium coating
- Anti-Retraction Valve

RRP	\$468
SALE PRICE	\$234
SAVE	\$234

Non-optic options available. Other connections offered for KaVo, Sirona, W&H, Bien Air.

All offers as listed. Valid until the 30th September 2023. Not to be used in conjunction with any other NSK promotions. All prices are as quoted in AUD and include GST.

## Surgical Handpieces

### Implant and Oral Surgical Handpieces

**OFFER 1**  
**BUY 1 GET 15% OFF RRP**

**OFFER 2**  
**BUY 2 GET 20% OFF RRP**

**OFFER 3**  
**BUY 3 GET 25% OFF RRP**

MODEL **Ti-Max X-SG20L**

ORDER CODE **C1003**

- Optic
- Reducing 20:1
- Titanium
- CA Burs

RRP	BUY 1
\$2,720	<b>\$2,312</b>
BUY 2	BUY 3
<b>\$2,176</b>	<b>\$2,040</b>



MODEL **S-Max SG20**

ORDER CODE **C1010**

- Non-Optic
- Reducing 20:1
- Stainless Steel
- CA Burs

RRP	BUY 1
\$1,657	<b>\$1,408</b>
BUY 2	BUY 3
<b>\$1,326</b>	<b>\$1,243</b>



MODEL **Ti-Max Z-SG45L**

ORDER CODE **C1107**

- Optic
- Increasing 1:3
- Triple Spray
- 25mm FG Burs

RRP	BUY 1
\$2,221	<b>\$1,888</b>
BUY 2	BUY 3
<b>\$1,777</b>	<b>\$1,666</b>



MODEL **Ti-Max X-SG65L**

ORDER CODE **H1009**

- Optic
- Direct Drive 1:1
- HP Burs

RRP	BUY 1
\$2,371	<b>\$2,015</b>
BUY 2	BUY 3
<b>\$1,897</b>	<b>\$1,778</b>



MODEL **SGA-ES**

ORDER CODE **H263**

- Non-Optic
- Direct Drive 1:1
- 20° angulation
- 65mm HP burs

RRP	BUY 1
\$2,107	<b>\$1,791</b>
BUY 2	BUY 3
<b>\$1,686</b>	<b>\$1,580</b>



MODEL **SGS-ES**

ORDER CODE **H264**

- Non-Optic
- Direct Drive 1:1
- 65mm HP Burs

RRP	BUY 1
\$1,520	<b>\$1,292</b>
BUY 2	BUY 3
<b>\$1,216</b>	<b>\$1,140</b>



MODEL **SGA-E2S**

ORDER CODE **H265**

- Non-Optic
- Increasing 1:2
- 65mm HP Burs

RRP	BUY 1
\$2,472	<b>\$2,101</b>
BUY 2	BUY 3
<b>\$1,978</b>	<b>\$1,854</b>



MODEL **SGS-E2S**

ORDER CODE **H266**

- Non-Optic
- Increasing 1:2
- 65mm HP Burs

RRP	BUY 1
\$1,782	<b>\$1,515</b>
BUY 2	BUY 3
<b>\$1,426</b>	<b>\$1,337</b>



**OSSEO 100 20% OFF**

## Osseo 100

### Osseointegration measuring device

#### Removes doubt

- Osseo 100 measures implant stability and osseointegration to enhance decisions about when to load the implant. Especially important when using protocols with shorter treatment time and treating higher risk patients.
- By mounting a MultiPeg™ the measurement is made in a second. A numerical result appears on an LED screen and that indicates how well the implant has integrated.



MODEL **Osseo 100**

ORDER CODE **Y1004175**

RRP **\$3,267**

SALE PRICE **\$2,614**

SAVE **\$653**

MultiPegs sold separately –  
for full list scan the  
QR code.



## VarioSurg3

Perform Safer Bone Surgery with 20% off NSK VarioSurg3

### Safer Bone Surgery

VarioSurg3 Ultrasonic bone surgery system provides safe and precise cutting of bone, without damaging soft tissue. A variety of tips makes procedures easier and stress-free to help provide more precise and efficient treatment.



- Control Unit
- LED Handpiece with 2m cord
- FC-78 Foot Control
- Sterilisation Cassette
- Handpiece Stand
- Irrigation tube (5 pcs.) and other accessories
- Basic H-S Kit (H-SG1, SG3, SG5 SG6D, SG7D, SG11 and Tip holder)

### VARIO SURG3 20% OFF

MODEL	<b>VarioSurg3 complete set</b>
ORDER CODE	<b>Y1002726</b>
RRP	<b>\$12,468</b>
SALE PRICE	<b>\$9,974</b>
SAVE	<b>\$2,494</b>
MODEL	<b>VarioSurg3 complete set without foot control</b>
ORDER CODE	<b>Y1002248</b>
RRP	<b>\$10,060</b>
SALE PRICE	<b>\$8,048</b>
SAVE	<b>\$2,012</b>

ALREADY A SURGIC PRO USER?  
PURCHASE VARIOSURG WITHOUT  
FOOT CONTROL AND GET A  
**FREE\* LINK SET**

MODEL	<b>VarioSurg3 SG Link set (Stand &amp; Cable)</b>
ORDER CODE	<b>Y1002729</b>
RRP	<b>\$1,116</b>
SALE PRICE	<b>FREE*</b>
SAVE	<b>\$1,116</b>

Create Synergy in Implantology by linking the Surgic Pro and Variosurg together with SG Link Set



STORE YOUR UNITS SAFELY WITH  
NSK iCart Duo  
**20% OFF**


## iCart Duo

The NSK iCart Duo allows you to install both the NSK SurgicPro and Variosurg Systems as well as additional accessories on to the one easily transportable cart.

MODEL	<b>iCartDuo</b>
ORDER CODE	<b>S9090</b>

- Dimensions: H101.65cm
- Weight: 16.5kg
- Adjustable Height
- Lockable Wheels

RRP	<b>\$2,355</b>
SALE PRICE	<b>\$1,884</b>
SAVE	<b>\$471</b>




## Variosurg3 Tip Kits

GET STARTED WITH  
**15% OFF**



MODEL	<b>Sinus Lift Kit</b>	MODEL	<b>Basic H-S Kit</b>
ORDER CODE	<b>Y900689</b>	ORDER CODE	<b>Y1002775</b>
• Tip Holder		• Tip Holder	
• SG1, SG3, SG6D, SG9, SG10, SG11		• H-SG1, SG3, SG5, SG6D, SG7D, SG11	
RRP	<b>\$2,098</b>	RRP	<b>\$2,183</b>
SALE PRICE	<b>\$1,783</b>	SALE PRICE	<b>\$1,856</b>
SAVE	<b>\$315</b>	SAVE	<b>\$327</b>



SURGIC PRO & VARIOSURG3  
CABLES **15% OFF**

Increase your workflow with additional Surgic Pro & VarioSurg3 cables

MODEL	<b>VS3-LED-HPSC</b>
ORDER CODE	<b>E1133</b>
• Optic scaler with 2m cord for VarioSurg3	
RRP	<b>\$5,517</b>
SALE PRICE	<b>\$4,589</b>
SAVE	<b>\$828</b>



MODEL	<b>SGL70M</b>
ORDER CODE	<b>E1023001</b>
• Optic micromotor with 2m cord for Surgic Pro	
RRP	<b>\$4,464</b>
SALE PRICE	<b>\$3,794</b>
SAVE	<b>\$670</b>

MODEL	<b>SG70M</b>
ORDER CODE	<b>E1025001</b>
• Non-Optic micromotor with 2m cord for Surgic Pro	
RRP	<b>\$3,548</b>
SALE PRICE	<b>\$3,016</b>
SAVE	<b>\$532</b>



## Perio-Mate Biofilm Eraser

Perio-Mate, the best biofilm eraser for gentle periodontal treatment.

Extensive range of applications ranging from supragingival to subgingival use

By removing the Perio-Mate nozzle tip, the subgingival area up to 3 mm from the gum line that requires the most frequent maintenance can be treated. In addition, there is a powder flow rate adjustment function to match operational site conditions, to enable a wide range of possible treatments from fast precise removal of biofilm to gentle and efficient biofilm removal around implants and in the delicate periodontal pocket.



## Powder Pro Care, above and below the gums


Efficient professional care that's kind to teeth.

It is easy to install powder therapy alongside existing facilities.

- The 'handy type' can be connected to the chair unit during use in one action.
- The nozzle tip is designed for safety and operability.
- The powder flow channel is short and direct with low risk of blockage.
- Fast and effective biofilm removal.

### Powder flow volume adjustment

The Perio-Mate's convenient Powder Control Ring makes fine-tuning powder flow simple and intuitive with quick fingertip adjustments. In addition to easily controlling and stopping the flow of powder in this way, you can cleanse the treatment area from biofilm uninterrupted.




### Smooth 360-degree rotation for safe and accurate operation

The handpiece can easily be rotated using the fingertips even while performing complex operations. The nozzle tip allows easy access to difficult to reach areas, and does not adversely affect neighbouring teeth and contours.



### Slim nozzle facilitates operation in the narrow confines of the oral cavity

The long, narrow nozzle shape is designed to allow easy access to areas that require inserting the handpiece at an angle, such as between molars.



### Product maintenance

The user can easily disassemble the nozzle, handpiece and powder case. Powder clogging is unlikely when handled and maintained correctly but is easy to manage by following the maintenance guidelines and the relevant accessories included with the product. All the parts of the handpiece (other than the nozzle tip) can be repeatedly cleaned ultrasonically or sterilized in an autoclave.



## PERIO-MATE 20% OFF

MODEL **PER-PTL-P for NSK Coupling**

ORDER CODE **Y1002654**

- Treatment of periodontal pockets and peri-implantitis
- Eradication of biofilm in periodontal pockets
- Removal of plaque and discolouration
- NSK's innovative nozzle tip
- Precise, flawless powder flow
- Slim nozzle design
- Smooth 360° rotation
- Powder Control Ring

RRP **\$1,696**

SALE PRICE **\$1,357**

SAVE **\$339**

MODEL **PER-WH-P for W&H® Roto Quick**

ORDER CODE **Y1002657**

RRP **\$1,809**

SALE PRICE **\$1,447**

SAVE **\$362**

MODEL **PER-SR-P for Sirona® Coupling**

ORDER CODE **Y1002656**

RRP **\$1,809**

SALE PRICE **\$1,447**

SAVE **\$362**

MODEL **PER-KV-P for KaVo® Coupling**

ORDER CODE **Y1002655**

RRP **\$1,809**

SALE PRICE **\$1,447**

SAVE **\$362**

**BUY 1 GET 1 FREE**

## Perio-Mate

### Biofilm Eraser Powder

The 'Perio-Mate Powder' used is highly water soluble to prevent sub gingival retention and has a lower Mohs hardness than dentine. Glycine ranks 2 on the Mohs hardness scale, while dentine is 2 to 2.5. Average particle size 25 microns.

MODEL **Perio-Mate powder (2 x 160g Bottles)**  
ORDER CODE **Y900938G**  
RRP **\$313**  
SALE PRICE **\$313**  
SAVE **\$313**



## Simple and precise root canal treatment

ENDO-MATE TC2 & AT **20% OFF**

### Endodontic Motors

#### ENDO-MATE TC2

MODEL **ENDO-MATE TC2 (230V)**  
MP-F16R head

ORDER CODE **Y1001027**

- Large LCD display with 5 programs for different file systems
- Supports most major brands of Ni-Ti files
- Torque control facility ranging from 0.1 to 4.0 NCM depending on head:gear ratio

RRP **\$2,263**

SALE PRICE **\$1,924**

SAVE **\$339**

#### ENDO-MATE AT

MODEL **ENDO-MATE AT (230V)**

ORDER CODE **Y1002190**

- Fits all major brands of Files
- Torque control, with integrated auto-reverse function
- Compact stylish control/monitor pad for simple and user friendly operation

RRP **\$1,845**

SALE PRICE **\$1,568**

SAVE **\$277**



**BUY ENDO-MATE TC2 OR ENDO-MATE AT AND GET 50% OFF IPEXII**

### Apex Locator

#### iPEX II

MODEL **iPex II Complete set**

ORDER CODE **Y1002208**

- Accurate measurement of any root canal
- SmartLogic technology
- Operates in all conditions, including dry, wet and bloody canals

RRP **\$1,696**

SALE PRICE **\$848**

SAVE **\$848**



"The NSK Endomate TC2 has been my motor of choice for the past 10 years. It is a well balanced, cordless motor that delivers consistent torque and speed for the endodontic treatment."

**Dr Morgan Wu**  
Specialist Endodontist  
Bds (Otago), Dcd (Endo)(Melb)



### iCLAVE PLUS OFFER

**Buy iClave Plus and get 50% OFF water purity kit + FREE WI-FI data logger**

#### B-TYPE VACUUM AUTOCLAVE\*

MODEL **iClave plus**

ORDER CODE **Y1003077**

- 20% more sterilisation space than a conventional stainless steel chamber of the same size.
- Copper chamber minimises internal temperature fluctuations.
- Super fast cycle as quick as 18 minutes, including drying.
- One year warranty.
- Unit size: H428 x W445 x D532mm, chamber 18L.

RRP **\$10,692**



**PURITY WATER PURIFIER 50% OFF + FREE WI-FI DATA LOGGER**

MODEL **PURITY OR2**

ORDER CODE **Z1284001**

RRP **\$2,810**

SALE PRICE **\$1,405**

SAVE **\$1,405**



**NSK Oceania Pty Ltd** [www.australia.nsk-dental.com](http://www.australia.nsk-dental.com)

Product Information | Technical Support

Freephone: 1300 44 33 21 (Australia) Email: [info@nskcoceania.com.au](mailto:info@nskcoceania.com.au)